Datum revize / Dátum revízie / Date of revision: 19.12.2024, ver.02

CZ

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POPIS VÝROBKU

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Medical brush

DESCRIPTION

Medical brush is used by hospital staff to scrub their hands, fingers, and forearms prior to surgery or other intervention where a high degree of personal hygiene is required. The product may also be used to mechanically clean the intact skin. This product has a synthetic sponge surface and a soft bristle surface. It may also include a separate nail cleaner and/or a disinfecting substance. Medical brush is made up of two parts, polyethylene and sponge.

INTENDED USE

Used to clean hands, fingers and forearms before surgery or other procedures requiring a high degree of personal hygiene; may also be used to mechanically clean intact skin.

USAGE

- Open the pouch.
- Take out the scrub brush, put the disinfecting solution on the sponge and scrub the hands and arms thoroughly with the brush part for about 3 minutes.
- Clean the fingers nails with the nail file.
- Wash hands with clean water at last and dry the hands before or after surgery.
- Discard after single use

CONTRAINDICATION

There are none known contraindications

Datum revize / Dátum revízie / Date of revision: 19.12.2024, ver.02

CAUTIONS / WARNINGS

- The product use should comply with the requirements of the relevant operation specifications of the medical authority and laws and regulations.
- The product can only be used by qualified medical personnel.
- Don't use after the expiry date
- If the sterile package is damaged, the product might be damaged or contaminated, in which case its sterility cannot be guaranteed. Therefore, in case of damaged inner packaging, foreign body inside the packaging, do not use.
- The product is for single use, discard after use.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

The storage area for this product must be cool, dry, well-ventilated and clean.

SHELF LIFE

The shelf life of this product is 5 years.

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