CZ

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POPIS VÝROBKU

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Disposable sterile syringe

DESCRIPTION

The Sterile Syringe consists of a graduated hollow barrel and a movable piston with a plunger. One end of the barrel has a nozzle, which permits attachment to a hub and a needle. The nozzle has a centric configuration. The barrel and piston are made from polypropylene (PP). The plunger is produced from rubber. The product is sterilized with Ethylene oxide. The product is supplied sterile, non-toxic, pyrogen-free and latex-free.

SPECIFICATION

1 ml, 2 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml

INTENDED USE

Syringe operated by hand is used for clinical liquid suction, and injection after being charged with liquid.

USAGE

- Do not use in case of damaged inner packaging, detached protective jacket or foreign body.
- Tear open the package at the tear of the single package, take out the product.
- In case of use with needle, attach the needle, remove the needle cap.
- According to conventional methods, pull the piston, suck in medical fluid, and immediately use.
- When withdraw blood or injection, use with the hypodermic needle, clamp the needle holder with a medical forceps, adjust the angle of the cutting edge, and rotate it 90° with a suitable force.
- Carry out routine operation; extract the patient's blood or inject

CONTRAINDICATION

- The product contains synthetic rubber, do not use if allergic to synthetic rubber.
- The product is prohibited from High-pressure injection, and the drugs that indicate incompatible with polypropylene and isoprene rubber hypodermic syringe.
- The product cannot be used with chloral and other incompatible chemicals

CAUTIONS / WARNINGS

- The product use should comply with the requirements of the relevant operation specifications of the medical authority and laws and regulations.
- Product can only be used by trained medical staff or under the guidance of medical staff.
- The product can only be used by qualified medical personnel.
- Don't use after the expiry date.
- If the sterile package is damaged, the product might be damaged or contaminated, in which case its sterility cannot be guaranteed and might even pose potential safety hazard.
- Therefore, in case of damaged inner packaging, foreign body inside the packaging, do not use.
- When the product is used with injection needle, screw the conical joint between the two tightly to prevent falling off or leakage.
- The product is for single use, discard after use.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

- Products in the transport process should prevent heavy pressure, direct sunlight, rain and snow dip.
- This product should be stored in the relative humidity of not more than 80%, non-corrosive gases, cool, dry and well-ventilated clean indoor.

SHELF LIFE

The shelf life of this product is 5 years.

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