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PRF20 – Paraffin wax, 20 kg

Intended use:

Used for infiltration and/or embedding of tissues or clinical specimens.

Intended user:

Healthcare professionals, laboratory professionals

Paraffin wax from MEDILAB ČR s.r.o. is intended for professional use only for in vitro diagnostics.

Description.

Highly pure and high quality paraffin with a melting point of 56-58 °C, with the addition of synthetic polymers, without DMSO (dimethyl sulfoxide) additives

Convenient pearl shape in white/colourless colour for a quick and easy installation process

Ensures smooth cutting, minimises compression and cell disruption

Cutting thickness reduced to 2µm.

Principle of the test:

Paraffin wax is used as an embedding medium to help preserve and prepare tissue samples for microscopic examination/microtomy.

Description of calibrators and controls, reagents:

Paraffin wax from MEDILAB ČR s.r.o. does not require the use of calibrators or reagents.

Sterility:

Paraffin wax is not sterile.

Instructions for use:

Wear protective gloves when using paraffin wax.

Paraffin wax should be used according to the guidelines in place at the medical facility/laboratory.

Paraffin wax is ready for use as soon as it is completely melted.

During processing, specimens should be completely saturated with paraffin wax for best results.

Analytical Functionality

Paraffin wax is not used for the detection of any particular analyte or marker. Paraffin wax is used to infiltration and/or embedding samples for microtomy. The functionality of paraffin wax from MEDILAB ČR s.r.o. analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitative evaluation, measurement range, linearity, cut-off, including the determination of appropriate sampling and handling criteria and the management of known relevant endogenous and exogenous interferences and cross-reactions, are not affected.

Clinical performance:

Paraffin is not intended to detect specific diseases or pathological processes and conditions.

Indicators of clinical performance such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, probability, as well as expected values in normal and affected populations do not affect the use of paraffin in the clinical setting.

Precautions

For use in in vitro diagnostics only, read the Material Safety Data Sheet for hazards and safe handling practices before use.

May only be used by appropriately trained and qualified personnel.

Before using the product, check for expiration date and do not use products with an expired expiration date.

For your health and safety, use this product in a well-ventilated area and use protective gloves and other safety precautions.

For single use only, do not reuse.

Report:

All serious adverse events occurring in connection with the product concerned must be reported to MEDILAB ČR s.r.o., as the manufacturer, and to the competent authority of the Member State in which the user is established.

Storage conditions:

Store in a clean, dry and well-ventilated place at a temperature up to 38°C. Protect from light.

Disposal:

Paraffin wax from MEDILAB ČR s.r.o. does not contain any infectious material. However, the samples must be treated as if they could transmit infection before and after fixation, and all materials that have come into contact with them must be disposed of according to the relevant safety measures for the healthcare facility.

Paraffin wax should be disposed of in accordance with the regulations of the healthcare facility/laboratory.

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