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Instruction for use**EN****DESCRIPTION**

Latex Foley Catheter is inserted through the urethra during urinary catheterization and into the bladder to drain urine from the bladder. It consists of shaft, drainage funnel, inflation funnel, balloon and valve.

The product can be used on perioperative patient, patient have urinary retention or bladder outlet obstruction, patient who need to measure the urine volume accurately.

FEATURES

- Made from latex which is soft and comfortable
- Silicone coated
- Col or coding for size identification
- Available with different types

CONTRAINDICATIONS

- Acute urethritis
- Acute prostatitis
- Acute epididymitis
- Allergy to latex

INTENDED USE/INDICATIONS

Latex Foley Catheter is inserted into the bladder through urethra and are indicated for routine urine drainage or for routine post-operative drainage and irrigation bladder.

ADVERSE REACTIONS

- Bleeding
- Urethral swelling
- Urethritis
- Catheter-associated urinary tract infection (CAUTI)

DESCRIPTION FOR USE**Preparation**

- The following additional items might be required for catheterization.
 - Sterile field, Sterile gloves, Items required for cleaning the patient meatus, (sterile, based on established techniques), Luer syringe with sterile water for balloon inflation, empty syringe for balloon deflate, sterile dressings.
 - Urine drainage device.
- Before use, visually inspect the catheter to detect mechanical damages and check that no leakage occurs.
- Place male patients in supine position, female patients in lithotomy position.
- Wash and dry hands thoroughly.

- Using aseptic technique, remove the catheter from its pouch and place it on a sterile field.
- Put on sterile gloves and remove the catheter sleeve.
- Lubricate the tip and the shaft of the catheter with a water-based lubricant.
- Clean the opening of the urethra and the surrounding area, using established techniques.

Catheterization

- With the uncontaminated hand, pick up the catheter from the sterile field.
- Carefully insert catheter tip into the bladder (normally indicated by urine flow), and then a further 3 cm to ensure balloon is also inside it.
- For catheters with Tiemann tip, the tip should be facing towards the ceiling or patients' face, assuming the patient is supine.
- Advance the catheter further to ensure that the balloon is located beyond the bladder neck, within the bladder.
- Inflate the balloon with the required fluid volume (marked on the catheter funnel).
- Slowly retract the catheter until some resistance is felt to ensure that the balloon is correctly located within the bladder, before the bladder neck.
- Connect the catheter to a urine drainage device.
- Observe urine flow.

Catheter removal

- For removal of the catheter, deflate the balloon by inserting a Luer syringe into the inflation valve. Release the syringe plunger and allow the balloon to deflate. Only use gentle aspiration, if necessary, to deflate the balloon. The balloon should spontaneously deflate.
- Discard the catheter according to hospital protocol.

Catheter care:

- Ensure local cleaning and hygiene protocols are followed to keep the catheter and meat us as clean as possible.
- Conduct hand hygiene immediately before and after any manipulation of the catheter and the urine drainage system. Wear disposable gloves when handling the system.
- Maintain unobstructed urine flow.
- Ensure that the catheter and collection tube do not kink.

Warnings/Precautions

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.

- The product is intended for use by physicians trained and aseptic technique must be practiced.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single use only.
- Not compatible with MRI for plastic valve catheter.
- This product contain Latex and might cause allergy. Doctors need to decide if the patient can use latex catheters according to patient's allergic history.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Do not over inflate the balloon. Refer to outer unit pack or funnel of catheter for balloon capacity.
- Inflate catheter balloon with sterile distilled water.
- Do not use ointment and lubricants having petroleum base. Use silicone oil or water soluble lubricants.
- Do not clamp catheter shaft. It may damage catheter and cause non deflation.
- Do not use needle syringe to inflate balloon. Use Luer syringe without needle.
- The retention time of Latex Foley Catheter should be decided by the medical care professionals according to the actual situation of using. Suggest to use less than 14 days.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

Storage conditions

- Store in a dark, cool and dry place.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

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